



MediaDent® Mercury Messenger Installation & Configuration

Server Installation

- From the server, download MediaDent® Mercury Messenger from the MediaDent® web site and save to your desktop.
- Once the download is complete, double click ServerSetup.exe to begin the installation.
- Click Next on the MediaDent® Mercury Messenger Server screen.
- Click yes if you accept the license agreement.
- By default MediaDent® Mercury Messenger is installed to C:\Program Files\MultiMedia Dental\MDMM, but this can be changed by clicking the Browse button and navigating to the desired location.
- Once you have reached the Start Install screen click Next to install MediaDent® Mercury Messenger.
- When prompted, click Finish to complete the installation.

Registration

- Once the installation is complete, you will be prompted to enter a registration key.
- To obtain a registration key, send an email to keygen@mediadentusa.com with MediaDent Mercury Messenger in the subject line.
- You will then receive a reply containing your MediaDent® Mercury Messenger activation key.
- Once you have received your activation key, click the Enter Key button.
- Enter the activation key provided into the appropriate fields and click OK.
- Once the Activation key has been entered you will receive a message that the Key is valid, and has been stored.

Important: Do not discard the activation key, as you will be required to enter this information on each workstation.

Server Configuration

- Once the activation key has been successfully entered, the MediaDent® Mercury Messenger configuration screen will appear.
- To configure or edit button captions, double click the desired button and enter the appropriate caption, then click Save.
- Once you have configured all button captions, click Save Buttons.
- You now need to add workstations. To do so, enter the first computer name, such as server, then click the add button. This will need to be done for each computer on the network.



- Once you have added all computers, you will then need to enter Personnel. To do so, type the staff members name and click Add. This will need to be done for each staff member.
- Once you have added all personnel, you will need to select a color for each staff member. To do so, click the drop down arrow located next to Edit Color and select the desired staff member.
- Once the staff member is active, click the Pick Color button, select the desired color, click OK, then click Save. Select the next staff member from the drop down menu and follow the same procedure. Do this for each staff member.
- If you have the need to delete workstations or personnel, use the drop down menu in the appropriate section, to select an item. Once the appropriate item is selected, click the delete button.
- Server configuration is now complete, click the exit button to close this window.

Sharing the MDMM directory for network access

- You will now need to share the MDMM directory. To do so, open windows explorer and navigate the location in which MediaDent® Mercury Messenger was installed (C:\Program files\MultiMedia Dental\MDMM by default)
- Right click on the MDMM folder and select Sharing and Security.
- Click the radio button next to Share this folder, and then click the permissions button.
- Allow Everyone full control by placing a check in the Allow box, next to Full Control, then click Apply, then OK.
- Click Apply on the MDMM properties window, then OK.
- Close windows explorer.

Workstation Installation

- From a workstation, open Windows Explorer and navigate through Network Places to the MDMM folder on the server.
- Locate WSsetup.exe and double click to begin the installation.
- Click Next on the Welcome Screen.
- You will then need to enter the server path. To do so, click the Browse button and navigate through network places to the MDMM folder. Highlight this folder and click OK.
- Click Next to begin the installation.
- When prompted, click the Finish button to complete the installation.

Registration

- Once the installation has completed, double click the MediaDent® Mercury Messenger Icon on your desktop.
- Upon launching, you will be prompted to enter the activation key.
- Using the Activation key received earlier, enter the information into the appropriate fields and click OK.



MultiMedia Dental Systems, Inc.
Engineering and Support Group
137 N. Broad St., Ste. E
Winder, GA 30680
800-631-5867 Fax: 770-867-5087
www.mediadentusa.com

- Once the Activation key has been entered you will receive a message that the Key is valid, and has been stored.

Workstation Configuration

- Once the key has been stored, you will be prompted with a configuration screen.
- Using the drop down menu located next to This Machine: select the computer name in which you are currently working on, then click Save.
- Workstation configuration is now complete. You can minimize MediaDent® Mercury Messenger to the taskbar, but do not close. If the application is closed, you will not receive messages.